



Quality Antibody Production Since 1960

Pocono Rabbit Farm & Laboratory Inc.

Dutch Hill Road, P.O. Box 240, Canadensis, PA 18325

Phone: (800) 622-6381 • Fax: (570) 595-9277

E-mail: antibody@prfal.com

Standard Guinea Pig Polyclonal Antibody Protocol for Fusion Protein Antigens

Subject	VAF Dunkin Hartley Female Guinea Pig
Adjuvant	Complete Freund's Adjuvant (CFA) initially, followed by Incomplete Freund's Adjuvant (IFA) for all subsequent injections.
Immunogen	Typically proteins and conjugated peptides. Typically 25ug to 200 µg immunogen is used for each immunization (however amounts may vary do to antigen availability).
Procedure	The immunogen is diluted to 0.125 ml with sterile saline and combined with 0.125 ml of the appropriate adjuvant. The antigen and adjuvant are mixed thoroughly to form a stable emulsion that is injected intramuscularly (IM) and/or intradermally into the hind haunches near the rear legs for a maximum of two times intradermally. This provides enhanced immune response from the sustained presentation of the immunogen. PRF&L has established immunization SOPs and Guidelines to be followed.
Bleeds	Blood is collected from the retro-orbital sinus by a microhematocrit tube into heparin. The plasma is decanted and clarified by centrifugation at 2,500 RPM for 15 minutes. PRF&L has established blood collection and laboratory SOPs and Guidelines to be followed.
Schedule	<p>Day 0 Preimmune Bleed yields about 3 to 5 ml of plasma. Immunize with 50 -- 100 µg/guinea pig antigen in CFA intradermally. (if concentration of the antigen is less than 0.8mg/ml also inject subcutaneous or IM)</p> <p>Day 28 Immunize with 25 --50 µg/guinea pig antigen in IFA intradermally. (if concentration of the antigen is less than 0.8mg/ml also inject subcutaneous or IM)</p> <p>Day 42 Immunize with 25 --50 µg/guinea pig antigen in IFA intramuscularly.</p> <p>Day 56 Bleed (3 – 5 ml plasma) Immunize with 25 -- 50 µg/guinea pig antigen in IFA intramuscularly Immunize with 25 -- 50 µg/guinea pig antigen intravenous (no adjuvant).</p> <p>Day 63 Bleed (3 – 5 ml plasma)</p> <p>Day 65 Ship preimmune and two bleed plasma to project sponsor for testing or PRF&L can perform ELISAs at the sponsor's request. The project can be discussed with PRF&L Staff and alterations to the schedule may be implemented at this point.</p> <p>Day 70 Immunize with 25 -- 50 µg/guinea pig antigen in IFA intramuscularly and continuing every 28 days.</p> <p>Day 84 Bleed (3 – 5 ml plasma) and, continuing every 28 days.</p>
Continuing	The guinea pig is maintained on a four-week cycle of immunizations and bleeds. The results of testing will determine whether the protocol should be continued, modified, or terminated. Extra weekly production bleeds and immunizations can be performed upon sponsor's request providing they fall within PRF&L Guidelines and SOPs.
Endpoint	The endpoint of the project is determined by the feedback received from the sponsor. He/she may request to terminate the project with serum (by exsanguination) or without serum (humanely euthanize or discard). PRF&L reserves the right to terminate a guinea pig prematurely that shows signs of morbidity.